



# Electric tables range

NIDEK offers a complete range of elegant and functional electric tables. Designed and manufactured in Saint-Priest, the Up tables are the outcome of NIDEK's expertise in the field of consultation furniture.

The variety of tabletop options meet all requirements and adapt to all spaces even the smaller ones.

All Ups offer a travel range of 66 cm (60 to 126 cm) and can be used while sitting or standing.



The new digital control panel indicates the tabletop height in real-time. Practitioners can set up to four preferential heights.



The table stand and castors guarantee a perfect mobility and keep balance. A single outlet (electrical and LAN plugs) is easily accessible.



## ■ Up mini

Perfect for smaller spaces, the Up Mini has a 1-instrument tabletop with compact footprint.



## ■ Up: 1 instrument table

The shape of this 1 instrument tabletop has been exclusively designed to fit between the armrests of wheelchairs. It comes as close as possible to disable patients, thus guaranteeing comfortable examination conditions.



## ■ Up: 2 instruments table

The Up 2-instrument table offers high stability and ergonomics. It is available in 2 tabletop shapes: V-top or straight top.



## ■ A wooden finishing touch

Practitioners can customize tabletops with a wooden touch. This brings unique and high-quality standards to their practice.




# Caractéristiques techniques\* Up

	Up Mini	Up: 1 instrument table	Up: 2 instruments V-shape	Up: 2 instruments straight-shape
References	HTA0018	HTA0019	HTA0020	HTA0021
Max. Size (with foot)	L 550xW 450 mm	L 952xW 563 mm	L 960xW 556 mm	L 1100xW 450 mm
Nb devices	1	1	2	2
Mini / Maxi height	60 cm / 126 cm			
Max loaded	80 Kg	80 Kg	2*40 Kg	2*40 Kg
Integrated connectivity	1 power supply / 1 LAN	3 power supply / 1 LAN		
Weight	34 Kg	44 Kg		
Power supply	100-230 VAC / 50-60 Hz			
Nominal power	1000 VA			


\*Specifications and design are subject to change without notice for improvement.

## ■ Dimensions


Up mini




Up: 1 instrument table



Up: 2 instruments V-shape



Up: 2 instruments straight-shape



**Indications:** Class I / CE medical device. The UP electric table (standard, PMR, 2 instruments V-tray, GlassPop) allows to perform an ophthalmologic examination. This device avoids the movement of the ophthalmologist and the patient, by adapting the position of the ophthalmology instruments to the latter electrically.

**Indication of use:** medical device intended for use by health professionals or by users trained by these health professionals. The safety precautions and procedures for use must be fully understood before using this device.

Please read the instructions in the user manual carefully.

Equipment manufactured by NIDEK SA. Last update date: september 2021.



**NIDEK SA**  
**Head-office**  
 Ecoparc  
 9, rue Benjamin Franklin  
 94370 Sucy-en-Brie - France  
 Tel.: +33 (0)1 49 80 97 97  
 Fax: +33 (0)1 49 80 32 08  
 Web: www.nidek.fr

**NIDEK SA**  
**Production facility**  
 Multiparc de Parilly  
 50, rue Jean Zay  
 69800 Saint-Priest - France